

Quantitative Reactions & Analysis

9/22/24

Isabella Guajardo

Elena Torres

Purpose: Learn how to set up and complete the filtration of a mixture, how to use volumetric glassware and pipettes, and how to use LoggerPro and a Spectrometer. Predict products from mixing two reactants, observe reaction, and analyze two methods used to determine the concentration of a solution. Determine which method is more accurate.

Procedure: Supplies/chemicals

- Water
- 0.1000 M CuSO_4
- 0.1000 M Na_2CO_3
- unknown M CuSO_4
- 50 mL beakers
- 100 mL beaker
- 250 mL beaker
- 5 mL volumetric pipet
- 10 mL volumetric flasks
- 10 mL graduated cylinders (in fume hood)
- pipet bulb
- watch glass
- funnel
- labeling tape
- filter paper
- 6 cuvettes
- SpectroVis Plus spectrometer

Part A: Reaction between Sodium Carbonate & Copper(II) Sulfate

- ① Get 10 mL of CuSO_4 (UNKNOWN) solution and 10 mL of the 0.1000 M of Na_2CO_3 solution using the labeled graduate cylinders in the fume hood.
- ② Get a 50 mL beaker and pour both solutions into it while slightly swirling it.
- ③ Get filter paper and record mass on analytical balance.
- ④ Fold filter into quarters and open it inside a funnel.